


How to assemble

There are two pieces that make the face shield: a transparency film and a 3D-printed “3DVerkstan” headband.

To assemble, align the holes in the transparency film with the protrusions on the headband, then push through. Band may be tray up or tray down, either is functional. Sometimes the film rips when attaching, so we’ve included extra sheets.



attach *transparency film*

to *headband*
(center holes first)

completed face shield

DISCLAIMER

These face shields are being donated free of charge to assist the recipients in dealing with the emergency situation created by the COVID-19 pandemic. They were manufactured in a home environment by non-specialists in the field without prior experience in the manufacture or design of such products, and without reference to any guidelines, instructions or treatises regarding the manufacture or design of face shields. They were also manufactured without reference to any federal, state or local regulations governing the production or design of such masks.

The face shields are being provided as is without any guarantees or warranty. In making this donation, we make no warranties of any kind, either express or implied, including but not limited to warranties of merchantability, fitness for a particular purpose, of title, of noninfringement of third party rights or of compliance with federal, state or local regulations. Use of the product by a user is at the user’s risk.

NO REPRESENTATIONS OR WARRANTIES, EITHER EXPRESS OR IMPLIED, OF MERCHANTABILITY, FITNESS FOR A SPECIFIC PURPOSE, NON-INFRINGEMENT OF INTELLECTUAL PROPERTY RIGHTS OF OTHERS, COMPLIANCE WITH FEDERAL, STATE OR LOCAL REGULATIONS OR OF ANY OTHER NATURE ARE MADE WITH RESPECT TO THESE PRODUCTS.

thank you for
your hard work
and bravery